


Bauerfeind GenuTrain® P3 Knee Support



GenuTrain P3 is a support that centers the patella during joint movement.

Product Features:

- Anatomically knitted support incorporating contoured silicone insert surrounding the patella
- An integrated guide keeps the patella from shifting out of place
- During activity the insert not only counteracts patellar malalignment but provides intermittent compression for the reduction of swelling and edema
- Promotes proprioception, thus heightening sensory awareness in the knee for increased joint stabilization
- Innovative trapezoidal stitch in the popliteal region prevents bunching and pinching
- Anatomical knit with graduated compression graduated toward the edges prevents constriction of circulation and rolling of the support
- Silicone banding at top of support and medial/lateral stays guard against migration
- Lightweight, breathable knit will not retain heat and is completely machine washable

Can be used to treat:

- Patellar subluxation
- Chondromalacia
- Patellar dislocation
- Osgood-Schlatter's disease after lateral release surgery
- Anterior knee pain
- Patellar instability
- Osteoarthritis

GenuTrain P3

Size	Loose Circumference in inches		Item #	
	6" Below knee	6 3/4" Above knee	(Right)	(Left)
1	11 - 12 ^{1/4}	15 ^{1/4} - 16 ^{1/2}	110414020*0601	110414020*0701
2	12 ^{1/4} - 13 ^{1/2}	16 ^{1/2} - 17 ^{3/4}	110414020*0602	110414020*0702
3	13 ^{1/2} - 14 ^{1/2}	17 ^{3/4} - 19	110414020*0603	110414020*0703
4	14 ^{1/2} - 15 ^{3/4}	19 - 20	110414020*0604	110414020*0704
5	15 ^{3/4} - 17	20 - 21 ^{1/4}	110141020*0605	110414020*0705
6	17 - 18 ^{1/4}	21 ^{1/4} - 22 ^{1/2}	110414020*0606	110414020*0706

*1=Nature, *8=Titanium

Point of measurement: Size is determined by loosely measuring:

Point 1. Circumference at 6 inches below center of the knee

Point 2. Circumference at 6^{3/4} inches above center of the knee

