

Bauerfeind MyoTrain® Thigh Support



MyoTrain is similar to a muscle taping bandage, but is faster, easier and provides improved results.

Directly after injury the additional acute pad can provide targeted compression to the injured area, to limit the consequences of the injury and shorten healing time. MyoTrain controls muscle tone and helps prevent repeated injury. It does not loosen with wear and pressure can be individually adjusted. MyoTrain does not slip, is washable, and can be re-used repeatedly.

MyoTrain offers innovative support for everyday or athletic use in the treatment of thigh injuries. It is a high-quality medical product that can be easily put on and removed.

MyoTrain works according to the principles of so-called muscle suspension bands (taping) using compression after muscle injuries and represents an effective therapeutic aid. The combination of anatomical, breathable knit, two adjustable compression pads . acute and functional pads . and individually adjustable straps provide the highest level of treatment and comfort.

Immediately following injury both compression pads are inserted into the support, to add targeted compression and to relieve injured muscles. The strap system provides individualized compression. The spread of the hematoma is minimized, and the extent of injury is limited and promotes the healing process. At a later stage, the MyoTrain can be worn without the compression pads.

Indications

- Muscle injuries to the thigh (especially cramps, pulls, ruptures, tears)
- Prevention, especially against repetitive injuries
- Bruising to thigh muscles
- Additional support for follow-up treatment of muscle injuries

Features

- The anatomical, breathable knit, with siliconized and pressure-reducing edges ensure an optimal fit and excellent comfort.
- Both pads (acute and functional pads) can be added and removed; MyoTrain can also be worn without pads
- The Velcro straps make it easy to put on and remove and provide individualized compression
- The anatomical knit also provides a proprioceptive effect, even when used without the straps.

Mode of action

- MyoTrain works according to the principles of so-called muscle suspension bands (taping) using compression.
- Immediately following injury, an additional, smaller acute pad can add targeted compression to stop hemorrhaging and the spread of hematoma
- The larger functional pad adds compression to the area around the injury to minimize the edematous swelling that typically arises.
- The lengthwise ribbing along the surface of the function pad improves lymphatic circulation.
- MyoTrain regulates muscle tone and reduces the risk of further injuries.

Size	1	2	3	4	5	6	7	Order number
Circumfer. ¾" above Knee	12 ½ - 13 ¾	13 ¾ - 15	15 . 16 ¼	16 ¼ - 17 ¼	17 ¼ - 18 ½	18 ½ - 19 ¾	19 ¾ - 21	1103250008000°
Circumfer. 1 1/4" above Knee	17 ¼ - 19	19 . 20 ½	20 ½ - 22	22 . 23 ½	23 ½ - 25	25 . 26 ½	26 ¾ - 28 ¼	

°=size