


Bauerfeind CaligaLoc® Stabilizing Ankle Brace



Use Bauerfeind CaligaLoc® to quickly and easily stabilize the entire ankle giving both medial and lateral support. CaligaLoc provides excellent ankle stabilization by counteracting supination (the outward rolling of the ankle) and preventing abnormal loads in the upper ankle. Can be worn around the clock, as a support splint in normal shoes during the day and as a positioning splint without shoes at night. Provides support of the foot in dorsiflexion. Keeps your ankle stable even if you have moderate foot drop and instability. Unlike stirrup supports, CaligaLoc stabilizes the ankle even when you are not wearing a shoe.

Product Features:

- Stabilizes the ankle by counteracting supination
- Prevents abnormal loads in the upper, and to a certain extent, lower ankle
- Anatomically contoured, thin, lightweight design ensures maximum comfort and convenience
- Heat moldable and can be reshaped if necessary for a individual customized fit
- Slightly pronating raised heel section (toward the outside) takes the strain off the lateral malleolar ligaments (outer ankle ligaments)
- Four adjustable Velcro straps provide easy application and proper fit
- Latex-free construction eliminates skin irritation
- Easy to clean with soap and warm water

Can be used to treat:

- Conservative treatment of torn ankle ligaments
- Charcot-Marie-Tooth Disease (CMT)
- Ankle stabilization
- Chronic ankle instability
- Drop-foot
- Post-operative protection after ligament suturing / reconstruction
- Post-Achilles rupture
- Night splint for plantar fasciitis

| Size | Shoe Size | | Item # | |
|------|-------------|------------|----------------|----------------|
| | Men | Women | Right | Left |
| 1 | 6 and below | 3 ½ - 7 | 12013003080601 | 12013003080701 |
| 2 | 6 ½ - 9 ½ | 7 ½ - 10 ½ | 12013003080602 | 12013003080702 |
| 3 | 10 - 14 ½ | 11 - 15 | 12013003080603 | 12013003080703 |